



**Laerdal**<sup>®</sup>  
helping save lives



## More Efficient CPR Training

### Resusci Anne<sup>®</sup> QCPR

The Resusci Anne<sup>®</sup> QCPR manikin and technology enables healthcare providers to obtain competence to deliver quality CPR in accordance with the latest Guidelines. QCPR technology is utilized to help both instructors and trainees more effectively monitor and review CPR performance, increasing the efficiency and value of training time.



[Laerdal.com](http://Laerdal.com)

# Resusci Anne® QCPR

In addition to providing high-quality CPR training, the Resusci Anne® QCPR manikin allows for those skills and others to be integrated into team training. The available add-on options and accessories make Resusci Anne QCPR more than just a high-quality CPR training manikin. The manikin can be upgraded to train teams on defibrillator use, airway management, IV administration, and rescue/trauma situations.

## Key Training Features

### CPR

- Quality CPR feedback (with feedback device)
- Mouth-to-mouth ventilations
- Ventilations with face shield, pocket mask, and BVM

### Airway

- Realistic chest rise/fall with ventilations
- Sellick maneuver
- Positive pressure ventilation
- Airway occlusion
- Supraglottic airway (with optional Airway Head)

### Pulses

- Manual Carotid pulses

### Defibrillation

- Train using a live defibrillator (with ShockLink™)
- Detect Pad Placement  
(with optional AED Skin, plus Laerdal AED Trainer or ShockLink)

### General

- Interchangeable pupils (with optional Airway Head)
- IV initiation and fluid administration (with optional IV Arm)
- Rescue/Trauma Simulation  
(with optional Trauma Limbs or Rescue Limbs)
- Latex free
- Wired and wireless feedback device options
- Bluetooth Low Energy
- Rechargeable battery (up to 34 hours operation time per charge)

## QCPR Feedback

- Compression rate and depth
- Complete release
- Correct hand position (not available on apps)
- Appropriate ventilation volume
- Amount of compressions and ventilations
- Frequency and length of interruptions (not available on SkillGuide/apps)


## Feedback Device Options



206-30001  
SimPad® PLUS with SkillReporter™  
206-50050  
SkillReporter Software License



170-30050  
SkillGuide™



202-56050  
SkillReporter (PC)



### Feedback Apps (free)

QCPR Learner  
   
QCPR Instructor

## Optional Limb Accessories



312025 /  
150-10150  
IV Arm



090050  
Bleeding  
Control Legs



312050  
First Aid/  
Trauma Limbs



312000  
Rescue Limbs

## Starting Configurations

### Resusci Anne® QCPR

171-00160	Torso
171-01260	Full Body
172-00160	w/ Airway Head – Torso
172-01260	w/ Airway Head – Full Body
173-00160	AED - Torso
173-01260	AED - Full Body
174-00160	AED w/ Airway Head - Torso
174-01260	AED w/ Airway Head - Full Body

### Kits

173-50750	AED Skin Upgrade Kit (add AED skin to non-AED skin configurations)
171-20050	Periodic Maintenance Kit
316000	Airway Trainer Update Kit (add Airway Head to non-Airway Head configurations)

### Warranty

One-year limited warranty. Refer to Laerdal Global Warranty for terms and conditions.

Information is subject to change. To learn more, visit [Laerdal.com](http://Laerdal.com)